

Finding out how much salt there is in seawater

In this experiment you will be using the microscale titration apparatus to find out how much salt there is in seawater. You will be titrating silver nitrate solution against seawater using potassium chromate as indicator.


Students must wear appropriate eye protection.

Instructions

1. Set up the microscale titration apparatus (see 'Apparatus and techniques for microscale chemistry' handout).
2. Fill the apparatus with silver nitrate solution (see 'Apparatus and techniques for microscale chemistry' and 'Element solutions' handouts).
3. Using a pipette add 1 cm³ of seawater to a 10 cm³ beaker.
4. Add one drop of potassium chromate indicator solution.
5. Titrate until a permanent red colour is observed and record the titre and any other observations.
6. Repeat and take the mean titre.
7. Calculate the salt content of the seawater. (The information below will help you with the calculations.)

Results and calculations

1. What do you think a reasonable answer for the percentage of salt in seawater is likely to be? Compare this with the value you get after working through the following calculations.
2. The equation for the reaction that occurs during the titration is:




What is the white precipitate that you see during the titration?

What causes the red colour at the end-point?

(Hint: think of some of the characteristics of transition element compounds)

3. From the equation it follows that:



Calculate the number of moles of silver nitrate used in the titration using the formula:

$$t/1000 \times C$$

where t = mean titre volume in cm³

and C = concentration of silver nitrate

4. Hence find the number of moles of chloride ion present in 1 cm³ of seawater.
5. Convert your answer in (3) to a mass by multiplying by the relative atomic mass of chlorine.
6. Convert your answer for (4) to a mass of sodium chloride and hence a percentage of sodium chloride in the seawater. What assumption are you making in your calculation?



7. Compare your result with your earlier guess. Are they similar?
8. Compare your results with others in your group.

Health and Safety

Students must wear appropriate eye protection (Splash resistant goggles to BS EN166 3).

Silver nitrate, 0.5 M AgNO_3 (aq) is CORROSIVE and will stain skin & clothing.

Potassium chromate, 0.2 M KCrO_4 is a carcinogen, mutagen and skin sensitiser as well as a skin/eye irritant. Explosive or vigorous-burning mixtures can be formed with aluminum and other metals and combustible materials.

Disposal: silver residues from a whole class should be collected and disposed of appropriately; small quantities can be washed to waste.

Credits

© Royal Society of Chemistry

Health & safety checked May 2018

Page last updated August 2018





© Royal Society of Chemistry, registered charity number 207890.